

Scientek™

Hospital and Laboratory Equipment

BOTTLE WASHING & FILLING

SW 1200A



SWBFS 25A

SW 1200 & SWBFS Series

SWBFS Bottle Filler



SWBFS 25A

Scientek Bottle Fillers provide uniform water delivery using either manual or automatic controls. Units are designed to accept common bottle basket configuration (24, 25 or 36 position). Custom sizes are also available. All Scientek Bottle Fillers are designed for ease of cleaning, low maintenance, and long life.

Standard Features Include:

- All stainless steel welded construction type 304 #4 polished finish.
- Easy glide white sanitary poly basket slides.
- Adjustable height manifold.
- Uniform volume no-drip nozzles.
- Swing away manifold access c/w quick water connection.
- Lower basket storage shelf.
- Auto models provide adjustable timed water fill and automatic reset c/w electronic solenoid valve operation.



SWBFS 2x36A



SW 1200 Bottle Washer

Scientek Bottle Washers have been designed to combine greater processing capacity with multi-cycle programming options. Choose from a single basket machine with automatic vertical load and unload doors, to a full indexing automated

multi-basket loading, washing and unloading module.

Whichever system suits your needs, they all come equipped with the reliable Allen Bradley Micro-Processor control, and Scientek's proven rotary spray washing system to deliver high impinging water pressure to get your products clean.

Standard Features Include:

- **Stainless steel conveyor belts.**
- **All stainless steel type 304 welded construction.**
- **Automatic vertical lift doors c/w viewing window.**
- **Multiple programmed cycle selection with user interface controls.**
- **Digital temperature control and guarantee.**
- **Powerful 7 1/2 HP stainless steel pump and motor.**
- **316 stainless steel steam sump heater.**



SW 1200A

SW1200 & SWBFS BOTTLE WASHER FILLER SERIES

Table of Scientek™ models to help you plan your facility.

MODEL #	INTERIOR/LOAD SIZE	EXTERIOR DIMENSIONS	BASKET SIZE
*SW 1200 Washer	19"W x 78"L x 15"H (483mm x 1981mm x 381mm)	34"W x 84"L x 83"H (864mm x 2134mm x 2108mm)	4 @ 19x19
*SW 1200 Auto Load	19"W x 76"L x 15"H (483mm x 1930mm x 381mm)	32"W x 78"L x 36"H (813mm x 1981mm x 914mm)	4 @ 19x19
*SW 1200 Auto Unload	19"W x 78"H x 15"D (483mm x 1981mm x 381mm)	32"W x 78"L x 36"H (813mm x 1981mm x 914mm)	4 @ 19x19
** SWBFS-25	25 Bottle Filler	22"W x 50"L x 48"H (559mm x 1270mm x 1219mm)	16x16
** SWBFS-36	36 Bottle Filler	25"W x 50"L x 48"H (635mm x 1270mm x 1219mm)	19x19

* Custom basket size configurations available to suit user needs and floor space.

** For automatic filler operation add suffix A to model number.

** For a change in bottle quantity add correct bottle number.

Bottle Washing Options

- Electric heating in place of steam
- Manual loading
- Automatic loading and unloading
- Number of baskets in load
- Basket size configuration
- DI or RO water final rinse
- Cycle data printing
- Chemical level monitoring
- PC reporting package
- On-line modem with factory
- Seismic anchoring

Bottle Filling Options

- Automatic timed filling
- 316 stainless steel construction
- DI or RO water supply
- Quick change manifolds
- Single hand filler
- Overhead bottle stopper station
- Dual filling manifolds
- Acid mixing station
- PH monitoring control
- Custom sizes & baskets
- Mobile unit with flexible utility connections

Scientek reserves the right to modify materials and specifications and to make product improvements without prior consent.

For additional product information please contact the factory or the Scientek representative in your area.

ScientekTM
Hospital and Laboratory Equipment

11151 Bridgeport Road
Richmond, BC, Canada V6X 1T3
Tel: (604) 273-9094
Fax: (604) 273-1262
E-mail: sales@scientek.net
www.scientek.net