



## BEEKAY RABBIT MAINTENANCE DIET (BK005P)

Beekay rabbit maintenance diet is suitable for rabbits maintained in laboratories for experimental purposes.

For breeding and young growing animals breeding diet is recommended.

### Product Information

Order Name: Beekay Rabbit Maintenance Diet (BK005P)

Main Ingredients: Wheat, Barley, Molasses, Dicalcium Phosphate, Grass Meal, , Extracted Linseed , Fat Blend, Minerals, Vitamins, amino acids

Form: Pellets, 3mm diameter

Presentation: 12.5kg paper sack.

Other packing available on request.

Irradiation: 25kGy, double wrapped and vacuum packed.

Other presentations on request.

### Labelling

Each bag is labelled with

- Product reference
- Expiry date
- Batch identification

### Quality Control Analysis

Independent laboratories perform physical-chemical and microbiological analyses on each batch of diet. Periodic full chemical analyses are conducted on selected batches to assure consistency over time.

### Calculated Composition

Crude Oil	%	3.88
Crude Protein	%	15.64
Crude Fibre	%	11.78
Ash	%	8.40
NFE	%	50.30
Dig Crude Oil	%	3.52
Dig Crude Protein	%	13.54
Gross Energy	MJ/kg	15.43
Dig Energy	MJ/kg	10.36
Vitamin A	iu/kg	11120.0
Vitamin D <sub>3</sub>	iu/kg	2102.0
Vitamin E	mg/kg	119.0
Thiamin	mg/kg	19.7
Riboflavin	mg/kg	24.8
Pyridoxine	mg/kg	20.2
Vitamin B <sub>12</sub>	µg/kg	51.5
Vitamin K	mg/kg	26.7
Folic Acid	mg/kg	11.5
Nicotinic Acid	mg/kg	82.8
Pantothenic Acid	mg/kg	34.5
Choline	mg/kg	1597.0
Inositol	mg/kg	2050.0
Biotin	µg/kg	473.0
Calcium	%	1.05
Total Phosphorus	%	0.55
Available Phosphorus	%	0.34
Magnesium	%	0.30
Sodium	%	0.31
Chloride	%	0.65
Potassium	%	0.97
Iron	mg/kg	210.0
Copper	mg/kg	17.0
Manganese	mg/kg	91.0
Zinc	mg/kg	91.0
Cobalt	µg/kg	580.0
Iodine	µg/kg	2090.0
Selenium	µg/kg	390.0
Lysine	%	0.77
Methionine	%	0.28
Linoleic Acid	%	1.13
Linolenic Acid	%	0.16